
















innovation begins here

providing specialized ophthalmic products & services



AOI
AMERICAN
OPTISURGICAL
INCORPORATED

TABLE OF CONTENTS

 ABOUT AMERICAN OPTISURGICAL INC.	
Message from the President	1
A Look at American Optisurgical Inc.	2
 PHACO SYSTEMS	
Vizual® Phacoemulsification System	3-5
Horizon Phacoemulsification Accessories	6
 PHACO HANDPIECES	
Compatible & New Phaco Handpieces.	7-8
 ACCESSORIES	
Phaco Tips	9-10
Administrative Tubing/Phaco Tubing Sets.	11-12
Phaco Accessories	13-14
I/A Handpieces & Accessories	15-16
Vitrectomy Cutters	17
Cautery Accessories	18
 INSTRUMENT CARE SYSTEMS	
QuickRinse™	19-21
 PHACO HANDPIECE REPAIR	
Repair & Refurbishment Services	22
 CONTRACT MANUFACTURING & GENERAL INFO	
Contract Manufacturing	23
General Information	24

MESSAGE FROM THE PRESIDENT



Herbert "Cam" Cameron
 President
 American Optisurgical Inc.



Dear AOI Customer,

This year we are celebrating an exciting 20-year landmark. We at American Optisurgical Inc. (AOI) are proud to provide specialized products and services to the ophthalmic community. AOI has grown to become a recognized leader in the ophthalmic industry by providing quality service and outstanding products to medical professionals around the world. We have our great customers to thank for our growth and success. On behalf of the entire staff at American Optisurgical Inc., I would like to express our personal appreciation to you.

As we move into the future, AOI will continue to provide professional, dependable service and the attention to detail you have come to expect from us. Our team of knowledgeable staff members keeps current on industry trends and developments in order to provide you with the best customer service and products available.

At AOI we look forward to working with you toward the ultimate goal of improving ophthalmic technology worldwide.

Herbert "Cam" Cameron

Herbert Cameron III

President, American Optisurgical Inc.

MISSION STATEMENT

AMERICAN OPTISURGICAL INC. IS DEDICATED TO THE RESEARCH AND TECHNICAL DEVELOPMENT OF OPHTHALMIC INSTRUMENTS AND INSTRUMENT CARE SYSTEMS WHILE PROVIDING SUPERIOR TECHNICAL SUPPORT AND SERVICE TO THE OPHTHALMIC COMMUNITY.





A LOOK AT AMERICAN OPTISURGICAL INC.

An ISO Certified Organization

Located in southern California, American Optisurgical Inc., also known as AOI, operates from a 10,000 square foot facility in the city of Lake Forest, California.

Herbert "Cam" Cameron founded AOI in 1992 for two reasons. First, he wanted to make a difference in the ophthalmic industry by providing customers with superior service at an affordable cost. Secondly, he wanted to achieve a personal dream of owning a company in an industry that helped people.

Under his direction, AOI has become a leader in the ophthalmic industry worldwide. Having achieved his initial goals, Cam's vision for the company has extended into designing and manufacturing innovative medical devices while continuing to grow AOI's business globally in providing superior service and products.

American
Optisurgical Inc.
26902 Vista Terrace
Lake Forest, CA 92630

Toll Free: 800.576.1266
Tel: 949.580.1266
Fax: 949.580.1270

Website: www.optisurgical.com

Email: info@optisurgical.com



Vizual[®] COMPACT POWER



The Vizual[®] Phacoemulsification System is a complete anterior segment surgical system. Lightweight and portable, it offers high-end performance and features in a small console.

Vizual[®] Features

- Small lightweight console
- Large foldable color touch screen
- User friendly controls
- Non-invasive vacuum sensing
- Multi-mode phaco
- Automated I/V pole interface
- USB upgradeable

- Ultra-light four crystal U/S handpiece
- Custom programs (up to 24 users)
- Easy loading non-invasive tubing cartridge

Vizual[®] Modes

- U/S, burst, micro
- Anterior Vitrectomy
- Coag (Diathermy)
- I/A



Vizual[®] Phaco System Standard Packages all include a one-year warranty.

551-0000-000 The Vizual[®] Phaco System Consists of: Phaco Unit, Foot Pedal, Handpiece and one of the packages shown below.

Package #1

551-7000-001 Travel Case, 30° AT SS Tip, Silicone Sleeves (2), Test Chamber, Collection Bag (10), Reusable Tubing & Cartridge Set, Fuse Set, Power Cord (**120V, US**), Operator's Manual, Vinyl Cover

Package #2

551-7000-002 All items in Package #1 plus Automated IV Pole Cart

Package #3

551-7000-003 Travel Case, 30° AT SS Tip, Silicone Sleeves (2), Test Chamber, Collection Bag (10), Reusable Tubing & Cartridge Set, Fuse Set, Power Cord (**220V, Europe**), Operator's Manual, Vinyl Cover

Package #4

551-7000-004 All items in Package #3 plus Automated IV Pole Cart



VIZUAL® PHACOEMULSIFICATION SYSTEM & ACCESSORIES

AOI Part #	Product Description
------------	---------------------

Components Sold Separately

553-0000-001	Vizual® Foot Pedal
300-9500-270	AOI 4-crystal Handpiece (Stainless with Titanium horn)
551-0000-001	Automated IV Pole Cart
430-2215-002	Fuse Set (set of 2)
430-4023-001	Power Cord (120V) USA
430-4024-001	Power Cord (220V) Europe
551-6001-001	Vizual® Operator's Manual
551-7001-001	Vinyl Cover for Vizual® unit
551-7002-001	Travel Case with Foam

Phaco Tubing Sets

551-0102-004	Reusable Tubing/Cartridge Set with 40 collection bags, non-sterile (box of 4)
---------------------	---

U/S Phaco Tips

American thread; 19 gauge; non-sterile; .036" ID / .042" OD; .9mm ID / 1.1mm OD
All tips below are smooth hub with tip wrench included

100-1000-001	0° round
100-1000-151	15° round
100-1000-301	30° round
100-1000-451	45° round
100-1000-601	60° round

American thread; 20 gauge (small bore); non-sterile; .029" ID / .035" OD; .74mm ID / .89mm OD

100-6001-001	0° round
100-6001-151	15° round
100-6001-301	30° round


American thread; 21 gauge (micro bore); non-sterile; .021 ID / .032" OD; .53mm ID / .81mm OD

100-6002-001	0° round
100-6002-151	15° round
100-6002-301	30° round

Phaco Accessories

All items are non-sterile except where noted.

600-1001-050	Administrative Tubing, sterile, disposable (pack of 50)
100-0050-285	American Silicone Sleeves, reusable (set of 10)
100-0050-292	American Silicone Sleeve, reusable, 20 gauge, small bore (set of 10)
475-1005-010	Test Chambers, reusable (set of 10)
475-1010-036	Mayo Covers, sterile (set of 36)



553-0000-001



300-9500-270



551-0102-004



430-4023-001



100-1000-301

VIZUAL® PHACOEMULSIFICATION SYSTEM & ACCESSORIES



551-7002-001



551-0000-001



551-0067-001

AOI Part #	Product Description
------------	---------------------

Phaco Accessories (cont.)

- | | |
|---------------------|--|
| 551-0067-001 | AOI Vizual® Collection Bags (set of 10) |
| 750-6019-001 | Wire Brush |
| 100-0025-018 | Tip Wrench, smooth hub, blue plastic, reusable (pack of 5) |
| 480-1017-001 | Autoclave Tray - Medium, with 3 inserts |

Irrigation/Aspiration Accessories

Titanium Handles & Tips

- | | |
|---------------------|--|
| 150-2020-540 | Titanium handle, standard, I/A round |
| 150-2010-513 | Titanium, 0 degree .3mm port |
| 150-2010-533 | Titanium, curved .3mm port |
| 150-2010-553 | Titanium, binkhorst, 150 degree .3mm port |
| 150-2010-623 | Titanium, 45 degree for silicone sleeve, .3mm port |

Stainless Handles & Tips

- | | |
|---------------------|---|
| 150-2020-510 | Stainless steel handle, standard, I/A round |
| 150-2010-113 | Stainless steel, 0 degree .3mm port |
| 150-2010-133 | Stainless steel, curved .3mm port |
| 150-2010-123 | Stainless steel, 45 degree .3mm port |
| 150-2010-143 | Stainless steel, 90 degree .3mm port |

NOTE: I/A tips with .4mm and .5mm ports are available by special order.

For additional Irrigation/Aspiration accessories, see pages 15-16

Vitrectomy Products

- | | |
|---------------------|---|
| 250-1700-026 | Reusable Vit Cutter, sterile, without sleeve (each) |
| 250-1700-030 | Infusion Sleeve, non-sterile, reusable (each) |

Cautery Accessories

- | | |
|---------------------|--------------------------------------|
| 475-1015-001 | Cautery Cable |
| 475-1035-001 | Forceps with angled tips |
| 475-1040-001 | Diathermy Jewelers Forcep |
| 475-1045-001 | Reusable Slim Slat Pencil (18 gauge) |
| 475-1046-001 | Reusable Slim Slat Pencil (20 gauge) |
| 475-1048-001 | Reusable Slim Slat Pencil (23 gauge) |

Miscellaneous Vizual® Phaco Accessories

- | | |
|---------------------|--|
| 551-6008-001 | Vizual® Service Manual |
| 553-1019-004 | Adhesive Foot Pedal Grip Set (Set of 4 strips) |





HORIZON PHACOEMULSIFICATION SYSTEM ACCESSORIES

AOI Part #	Product Description
------------	---------------------

Components Sold Separately

All items are non-sterile except where noted.

550-0050-001	Horizon Footswitch
300-9500-270	AOI 4-crystal Handpiece (Stainless with Titanium horn)
430-2215-002	Fuse Set (set of 2)
430-4023-001	Power Cord (120V) USA
430-4024-001	Power Cord (220V) Europe
550-6006-001	Horizon Operator's Manual
550-7001-001	Vinyl Cover for unit

Phaco Tubing Sets

550-0102-004	High-Efficiency, reusable with T-Fitting (box of 4)
---------------------	---

U/S Phaco Tips

American thread; 19 gauge; non-sterile; .036" ID / .042" OD; .9mm ID / 1.1mm OD
All tips below are smooth hub with tip wrench included.

100-1000-001	0° round
100-1000-151	15° round
100-1000-301	30° round
100-1000-451	45° round

For additional phaco tips, see page 9-10.

Phaco Accessories

600-1001-050	Administrative Tubing, sterile, disposable (pack of 50)
100-0050-285	American Silicone Sleeves, reusable (set of 10)
475-1005-010	Test Chambers, reusable (set of 10)
475-1010-036	Mayo Covers, sterile (set of 36)
475-1055-050	Collection Bags (set of 50)
100-0025-018	Tip Wrench, smooth hub, blue plastic, reusable (pack of 5)

For additional accessories, see pages 13-14.

I/A Handles & Tips - See pages 15-16

Vitrectomy Products

250-1700-026	Reusable Vit Cutter, sterile, without sleeve (each)
250-1700-030	Infusion Sleeve, non-sterile, reusable (each)

Cautery Accessories

475-1045-001	Reusable Slim Slat Pencil (18 gauge)
475-1040-001	Diathermy Jewelers Forcep
475-1035-001	Forceps with angled tips
475-1015-001	Cautery Cable

Miscellaneous Horizon Accessories

550-6007-001	Horizon Service Manual
---------------------	------------------------



Unit no longer available



550-0050-001



300-9500-270



550-0102-001



COMPATIBLE & NEW PHACO HANDPIECES

AOI Part

Product Description

When you choose an AOI Compatible U/S Phacoemulsification Handpiece you will reduce your operational costs significantly without compromising quality or performance. Manufactured under the ISO quality control system, AOI's compatible handpieces are designed to perform to required specifications. AOI guarantees each handpiece with a standard 1-year warranty.*

Compatible Phaco Handpieces

500-5000-100 Alcon Series Ten Thousand Compatible



300-5000-130 Alcon Series 20000 Legacy Compatible



300-5000-110 Alcon MacKool Compatible (reverse thread)



300-7000-100 Bausch & Lomb (Storz) Millennium Compatible



300-1000-100 Bausch & Lomb (Storz) Premiere Protégé



300-1200-100 Allergan Sovereign Compatible



300-1300-100 AMO Series 4 Compatible





AOI Part #	Product Description
AOI's Slim Compatible Handpieces	
300-8000-200	Oertli Compatible
300-9400-100	Mentor Mini
300-9500-240	40 KHz 4-crystal
300-9500-270	AOI 4-crystal (for use with Horizon/Vizual®)



AOI handpieces come standard with silicone over-mold strain reliefs at critical cable connection points to extend instrument life and performance. Every AOI Compatible Handpiece is thoroughly tested for quality and performance before leaving the facility.

Please note: If you are interested in a handpiece type that is not listed above, please contact an AOI Customer Representative for more information.

**Please contact AOI Customer Service to receive warranty information on compatible handpiece exchanges.*



PHACO TIPS

AOI Part #	Product Description
------------	---------------------

For use with: Bausch & Lomb Millennium, Premier, Protégé; AMO/OMS Phaco Plus, Profinesse, Diplomat, Diplomate, Diplomax, Opsys; Allergan Sovereign, Series 4; Mentor SIStem, Odyssey; Alcon Universal/10K; AOI Horizon, Vizual®; Oertli.

Details: American thread; 19 gauge; non-sterile; .036" ID / .042" OD; .9mm ID / 1.1mm OD
Tips below are smooth hub with tip wrench included.

- 100-1000-001** 0° round
- 100-1000-151** 15° round
- 100-1000-301** 30° round
- 100-1000-451** 45° round
- 100-1000-601** 60° round

For use with: Alcon Series 20000 Legacy, Accurus, Infiniti.

Details: American thread; 19 gauge; non-sterile; .036" ID / .042" OD; .9mm ID / 1.1mm OD
Tip wrench not included.

- 100-2000-001** 0° round
- 100-2000-151** 15° round
- 100-2000-301** 30° round
- 100-2000-451** 45° round
- 100-2000-601** 60° round

For use with: Bausch & Lomb Phacotron Gold; Optikon 2000 Phacovisc; United Surgical Systems Plus II.

Details: Metric thread; 19 gauge; non-sterile; .036" ID / .042" OD; .9mm ID / 1.1mm OD
Tips below are smooth hub with tip wrench included.

- 100-3000-001** 0° round
- 100-3000-151** 15° round
- 100-3000-301** 30° round
- 100-3000-451** 45° round
- 100-3000-601** 60° round

For use with: Optikon 2000 P4000, Pulsar, Quasar.

Details: Metric thread; 19 gauge; non-sterile; .036" ID / .042" OD; .9mm ID / 1.1mm OD
Tips below are smooth hub with tip wrench included.

- 100-3050-001** 0° round, .79mm ID, .97mm OD
- 100-3050-151** 15° round, .79mm ID, .97mm OD
- 100-3050-301** 30° round, .79mm ID, .97mm OD
- 100-3050-451** 45° round, .79mm ID, .97mm OD
- 100-3050-601** 60° round, .79mm ID, .97mm OD





AOI Part #	Product Description
100-3075-001	0° round, .91mm ID, 1.1 mm OD
100-3075-151	15° round, .91mm ID, 1.1 mm OD
100-3075-301	30° round, .91mm ID, 1.1 mm OD
100-3075-451	45° round, .91mm ID, 1.1 mm OD
100-3075-601	60° round, .91mm ID, 1.1 mm OD



For use with: Bausch & Lomb Millennium, Premier, Protégé; AMO/OMS Phaco Plus, Profinesse, Diplomat, Diplomate, Diplomax, Opsys; Allergan Sovereign, Series 4; Mentor SiStem, Odyssey; Alcon Universal/10K; AOI Horizon, Vizual®; Oertli.

Details: American thread; 20 gauge (small bore); non-sterile; .029" ID / .035" OD; .74mm ID / .89mm OD

Tips below are smooth hub with tip wrench included.

100-6001-001	0° round
100-6001-151	15° round
100-6001-301	30° round

For use with: Bausch & Lomb Millennium, Premier, Protégé; AMO/OMS Phaco Plus, Profinesse, Diplomat, Diplomate, Diplomax, Opsys; Allergan Sovereign, Series 4; Mentor SiStem, Odyssey; Alcon Universal/10K; AOI Horizon, Vizual®; Oertli.

Details: American thread; 21 gauge (micro bore); non-sterile; .021" ID / .032" OD; .53mm ID / .81mm OD

Tips below are smooth hub with tip wrench included.

100-6002-001	0° round
100-6002-151	15° round
100-6002-301	30° round

For use with: Bausch & Lomb Millennium, Premier, Protégé; AMO/OMS Phaco Plus, Profinesse, Diplomat, Diplomate, Diplomax, Opsys; Allergan Sovereign, Series 4; Mentor SiStem, Odyssey; Alcon Universal/10K; AOI Horizon, Vizual®; Oertli.

Details: American thread; FACT (Fluid Actuated Control Tip); non-sterile; .021" ID

Tips below are smooth hub with tip wrench included.

100-1006-301	30° 19 gauge, round
100-1007-301	30° 20 gauge, round

For use with: Bausch & Lomb Millennium, Premier, Protégé; AMO/OMS Phaco Plus, Profinesse, Diplomat, Diplomate, Diplomax, Opsys; Allergan Sovereign, Series 4; Mentor SiStem, Odyssey; Alcon Universal/10K; AOI Horizon, Vizual®; Oertli.

Details: American thread; 2021; non-sterile; .029/.021" ID, .035/.032" OD

100-1005-301	30° round, smooth-hub with wrench
--------------	-----------------------------------

*Elliptical tips are available by special order request only.
Not all phaco tips are carried in stock.*

PHACO TUBING SETS



600-1001-050



600-2001-012



600-3002-004



600-3005-001

AOI Part #	Product Description
------------	---------------------

Administrative Set

600-1001-050	5' 5" Length; Roller Clamp; 5 Micron Filter
---------------------	---

For use with: Bausch & Lomb Surgical Model: Phacotron Gold

600-2001-012	High efficiency, disposable, sterile (box of 12)
---------------------	--

600-2002-006	High efficiency, reusable, non-sterile (box of 6)
---------------------	---

For use with: AMO (formerly OMS) Models: Diplomat, Diplomate, Diplomax, Opsys

600-3002-004	Standard, reusable, sterile (box of 4)
---------------------	--

600-3003-012	High efficiency, disposable, sterile (box of 12)
---------------------	--

600-3004-004	High efficiency, reusable, sterile (box of 4)
---------------------	---

600-3005-001	T-fitting, reusable, non-sterile (each)
---------------------	---


For use with: Bausch & Lomb Surgical/Storz Models: Premier, Protégé, Millennium

600-5100-008	Basic - disposable, sterile, with cassette, administrative set (box of 12)
---------------------	--

600-5100-011	Disposable Cassette - sterile, with molded tube (box of 6)
---------------------	--



AOI Part #	Product Description
600-5000-004	Reusable Cassette - non-sterile, with molded tube (box of 4)
600-5001-004	Reusable, non-sterile, without cassette (box of 4)
For use with: AOI Horizon	
550-0102-004	High-Efficiency, reusable, non-sterile, with T-Fitting (box of 4)
For use with: AOI Vizual®	
551-0300-006	Disposable Tubing/Cartridge Set, sterile (box of 6)
551-0102-004	Reusable Tubing/Cartridge Set with 40 collection bags, non-sterile (box of 4)
For use with: Mentor Models: Odyssey, SiStem	
600-9004-002	SiStem reusable, non-sterile (box of 2 with 40 collection bags) (Old Mentor P/N 22-3021)
600-9002-002	Odyssey reusable, non-sterile (box of 2 with 40 collection bags) (Old Mentor P/N 22-1860)



600-5100-011



600-5001-004



551-0102-004



600-9004-002



PHACO ACCESSORIES



100-0025-012



100-0025-014



100-0025-018



100-0050-285



100-0050-292



425-4003-010



425-4004-010



475-1005-010

AOI Part #

Product Description

Autoclave Tray, Collection Bags, Covers, Luer Adapters, Phaco Tip Wrenches, Silicone Sleeves, Test Chambers, Tip Protectors

All items are non-sterile except where noted.

100-0025-012 Duo Tip Wrench, square hub/smooth hub, titanium, reusable (each)

100-0025-014 Tip Wrench, smooth hub, stainless steel, reusable (each)

100-0025-018 Tip Wrench, smooth hub, blue plastic, reusable (pack of 5)

100-0050-285 American Silicone Sleeves, reusable (set of 10)

100-0050-292 American Silicone Sleeve, reusable, small bore, 20 gauge (set of 10)

425-4003-010 Female Luer Adapters, reusable (set of 10)

425-4004-010 Male Luer Adapters, reusable (set of 10)

475-1005-010 Test Chambers, reusable (set of 10)





AOI Part #	Product Description
(cont.)	
475-1010-036	Mayo Covers, sterile (set of 36)
475-1055-050	Collection Bags (set of 50)
475-1065-010	OMS Tip Protectors, reusable (set of 10)
480-1017-001	Autoclave Tray - Medium, with 3 inserts
551-0067-001	AOI Vizual® Collection Bags (set of 10)
600-3005-001	OMS T-fitting, reusable (each)
810-1000-001	AOI CD Catalog (PDF format)




475-1010-036



475-1055-050



475-1065-010














480-1017-001



600-3005-001



I/A HANDPIECES & ACCESSORIES



AOI Part #

Product Description

Titanium Irrigation/Aspiration Handles and Tips

150-2020-540 Titanium handle, standard, I/A round

150-2010-513 Titanium, 0 degree .3mm port

150-2010-533 Titanium, curved .3mm port

150-2010-523 Titanium, 45 degree .3mm port

150-2010-543 Titanium, 90 degree .3mm port

150-2010-553 Titanium, binkhorst, 150 degree .3mm port

150-2010-613 Titanium, 0 degree for silicone sleeve, .3mm port

150-2010-580 Luer Adapter (irrigation)

150-2010-633 Titanium, curved for silicone sleeve, .3mm port











150-2010-623 Titanium, 45 degree for silicone sleeve, .3mm port

150-2010-590 Luer Adapter (aspiration)





AOI Part #	Product Description
Stainless Steel Irrigation/Aspiration Handles & Tips	
150-2020-510	Stainless steel handle, standard, I/A round
150-2010-113	Stainless steel, 0 degree .3mm port
150-2010-133	Stainless steel, curved .3mm port
150-2010-123	Stainless steel, 45 degree .3mm port
150-2010-143	Stainless steel, 90 degree .3mm port
150-2010-153	Stainless steel, binkhorst, .3mm port
150-2010-213	Stainless steel, 0 degree, for silicone sleeve, .3mm
150-2010-233	Stainless steel, curved, for silicone sleeve, .3mm port
150-2010-223	Stainless steel, 45 degree, for silicone sleeve, .3mm
150-2010-180	Luer Adapter (irrigation)
NOTE: I/A tips with .4mm and .5mm ports are available by special order. Prices and availability vary depending on style of tip.	
Irrigation/Aspiration Accessories	
150-2040-004	Push-on O-rings (for titanium I/A tips) (5 sets)
150-2040-002	Bayonette style O-rings (for SS I/A tips) (5 sets)
480-1017-001	Autoclave Tray - Medium, with 3 inserts



VITRECTOMY CUTTERS

AOI Part #	Product Description
------------	---------------------

Disposable - sterile, without sleeve

250-1700-040	Vit Cutter, 20 gauge, low pressure, anterior segment only (box of 6) (Horizon, Vizual®)
---------------------	---

250-1700-060	Vit Cutter, high speed, 20 gauge, posterior segment only (box of 6) (Mentor SIStem)
---------------------	---

250-1700-070	Vit Cutter, high speed, 23 gauge, posterior segment only (box of 6) (Mentor SIStem)
---------------------	---

Reusable - sterile, without sleeve

250-1700-011	Vit Cutter, 20 gauge, posterior segment only (each) (Mentor SIStem)
---------------------	---

250-1700-026	Vit Cutter, 20 gauge, anterior segment only (each) (Horizon, Vizual®)
---------------------	---

Vit Cutter Sleeve (for use with vit cutters shown above)


250-1700-030	Infusion Sleeve, non-sterile, reusable (each)
---------------------	---



250-1700-011



250-1700-030



Additional styles of vit cutters may be available upon request.



CAUTERY ACCESSORIES

AOI Part #

Product Description

All items are non-sterile.

475-1015-001 AOI Horizon, Vizual®, Optikon P4000 & Pulsar Cautery cable

475-1020-001 Standard Bipolar Cable

475-1025-001 OMS Cautery cable

475-1030-001 Bausch & Lomb Storz Premier Cautery cable

475-1031-001 Bausch & Lomb Storz Millennium (Erbe) Cautery cable

475-1032-001 Storz Daisy Cautery cable

475-1035-001 Forceps with angled tips

475-1040-001 Diathermy Jewelers Forcep

475-1045-001 Reusable Slim Slat Pencil (18 gauge)

475-1046-001 Reusable Slim Slat Pencil (20 gauge)

475-1047-001 Reusable Slim Slat Pencil (23 gauge)

475-1048-001 Reusable Bipolar Pencil, tapered tip (20/23 gauge) sharp

Please note: Additional cautery products are available upon request.

Please contact an AOI Customer Representative for more information.



475-1020-001




475-1025-001



475-1030-001



475-1040-001




475-1046-001



QuickRinse™ AUTOMATED INSTRUMENT RINSING SYSTEM

Easy and efficient at the same time.

The QuickRinse™ Automated Instrument Rinse System standardizes post-surgical care of your reusable instruments and tubing prior to sterilization.



(footpedal optional)

Patent No.
6,206,014 B1

- *Designed for rinsing phaco and I/A handpieces, cannulas, vit cutters, reusable tubing and any other micro-surgical instrument with lumen*
- *Replaces syringe rinsing methods*
- *Helps reduce damage to instruments from excessive handling*
- *Automated cycles for consistent rinsing and air drying with a manual override feature*
- *Compact and lightweight – unit has 8.75” X 8.75” footprint and weighs only 4 pounds*
- *Bottle, cap, and tubing disconnect easily for sterilization*
- **System includes** adapters, nozzles, bottle assembly and a cleaning brush
- **Optional accessories include:**
 - *A waterproof foot pedal offers “hands-free” operation*
 - *Nose cone flush attachment*
 - *Stellaris™ Irrigation line flush attachment*
 - *Dual rinse “Y” fitting attachment*
 - *Filtered air line*

QuickRinse™ Automated Instrument Rinse System

750-0000-000

QuickRinse™ unit, power supply, power cord, complete bottle assembly, accessory kit, wire brush, operator’s manual, and instructional DVD



QUICKRINSE™ SYSTEM & ACCESSORIES

AOI Part #	Product Description
750-2051-001	Footpedal
425-4003-010	Female Luer Adapters, reusable (pkg. of 10)
425-4004-010	Male Luer Adapters, reusable (pkg. of 10)
430-4023-001	Power Cord (120V) USA <i>(included in system)</i>
430-4024-001	Power Cord (220V) Europe
750-1100-001	Power Supply (Universal) <i>(included in system)</i>
750-3001-001	Bottle, polypropylene (1000ml)
750-3008-004	Bottle Cap Gasket (set of 4)
750-3010-001	Complete Bottle Assembly (Bottle, Cap, Tubing, and Fittings) <i>(included in system)</i>
750-3011-001	Bottle Cap Assembly (includes Tubing & Fittings)



750-2051-001



425-4003-010



430-4023-001



750-1100-001



750-3010-001




QUICKRINSE™ SYSTEM & ACCESSORIES



750-3014-001



750-4001-001



750-6019-001



750-7004-001



750-7000-001



750-7009-001

AOI Part

Product Description

750-3013-001	Air Supply Tubing Assembly (Green)
750-3014-001	Air/Fluid Tubing Assembly (Blue/Clear)
750-3019-001	Complete Tubing Set
750-4000-004	Instrument Interface Tubing (Set of 4)
750-4001-001	Agneta Dual Rinse "Y" Fitting (for rinsing of 2 HP's simultaneously)
750-6018-001	Operator's Manual <i>(included in system)</i>
750-6019-001	Wire Brush <i>(included in system)</i>
750-6023-001	QuickRinse™ Instructional DVD <i>(included in system)</i>
750-7000-001	Accessory Fitting Kit: Includes: (2) Female & (2) Male luer adaptors; (2) 16 gauge nozzles & (2) 22 gauge nozzles <i>(included in system)</i>
750-7000-002	Fluid/Air Tubing Coupling O-ring Set
750-7004-001	Nose Cone Flush Attachment
750-7009-001	Air Supply Tubing Assembly with Filter
760-0020-001	Bottle Cap (Fittings & Gasket only)
<i>The MK QuickRinse™ System is no longer available.</i>	
MKQR Accessories	
430-4023-001	Power Cord (120V) USA
430-4024-001	Power Cord (220V) Europe
760-0005-001	MKQR Bottle & Cap
760-0025-001	MKQR Nozzle (white)
760-7001-002	MKQR Microfilters (pkg. of 2)
760-7002-001	Operator's Manual
760-7003-025	MKQR Cleanser / Lubricant (box of 50 2cc ampules)

HANDPIECE REPAIR & REFURBISHMENT SERVICES

AOI provides 40 years of combined experience in the repair and refurbishment of phacoemulsification handpieces and units. Selecting AOI as your repair facility is a sound investment that will provide you with savings year after year. AOI is an ISO 13485 certified company, which reaffirms AOI's commitment to quality.

All repairs are done at our headquarters in Lake Forest, CA.

WE SERVICE PHACO HANDPIECES FROM MOST MAJOR MANUFACTURERS INCLUDING:

Alcon®	Dorc®	MST®
AMO (now Abbott Medical Optics)/OMS®	Fritz Ruck	Optikon 2000™
Bausch & Lomb Surgical/Storz®	Geuder®	Staar Surgical®
	Mentor	

Complete Handpiece Rebuild includes:

- New Ultrasonic Motor (crystals and electrodes)
- New Seals Throughout
- New Connector
- Strain Relief
- New Cable



Free Evaluation

An AOI technician will conduct a free, thorough performance evaluation of your instrument before starting any work. Only upon your approval will any work commence.

Standard Upgrades

Each refurbished phaco handpiece receives standard upgrades that will extend its life and performance. These upgrades are included in the price.

Reconditioned Equipment

If you are interested in purchasing Reconditioned Equipment, contact American Optisurgical Inc. for more information.

Other Repair Services

- Grieshaber Vit Cutters
- I/A Handpieces
- Vit Cutters
- Selected Phaco Units: OMS (all), Storz® Premier, Protege, Gold, Horizon & Vizual®

Repair/Exchange Program

Need a working handpiece right away? AOI offers a handpiece exchange program to get you back into surgery (based on availability).

Warranty Periods

AOI offers a variety of warranty periods to fit into your facility's budget. AOI offers cost saving customized service contracts with the flexibility to include any number of handpieces and/or units that your facility may have. These contracts have proven to be financially beneficial.

PRODUCT DEVELOPMENT STAGES

AOI offers product development and contract manufacturing services. Please call to arrange a consultation.




Custom design and assembly



Custom machined parts



Haas lathe and mill



Quality control



Controlled environment



Injection molding



Laser welding



GENERAL INFORMATION & POLICIES

COMPANY LOCATION

American Optisurgical Inc.
26902 Vista Terrace
Lake Forest, CA 92630

CONTACT US

Toll free	800.576.1266 (U.S. only)
Tel	949.580.1266
Fax	949.580.1270
Email	info@optisurgical.com
Website	www.optisurgical.com

BUSINESS HOURS

Monday through Friday
7:00am to 3:30pm (Pacific Time)
Closed weekends and most national holidays

GENERAL POLICIES

Pricing

All prices are subject to change without notice.

Payment

Payment for goods and services in USD only. AOI accepts cash, major credit cards, company checks and bank-to-bank wire transfers.

Terms

Domestic Terms: Payment terms are Net 30 days. An interest charge of 1.5% will apply to all balances over 30 days.

International Terms: All international orders require prepayment including all shipping charges (determined on a per-order basis).

Disclaimer

Alcon®, AMO® (now Abbott Medical Optics) /OMS, Bausch & Lomb Surgical / Storz®, Dorc®, Fritz Ruck, Geuder®, MST, Mentor®, Staar Surgical®, Optikon 2000™ may be registered names and trademarks of their respective companies and manufacturers. AOI does not imply representation or affiliation with any of these companies in any manner.

AOI has obtained registration to ISO:13485 standards. AOI also has many products that are CE marked.

All photos of products are not to scale.

AOI reserves the right to discontinue and/or update product.

RGA (Return Goods Authorization)

Any items sent for repair must be accompanied by an RGA. Product returns also require an RGA.

Returns

All products must be returned within 10 days of receipt of product. An RGA must accompany the product. Upon approval of return, a 20% restocking fee may be applied.

Warranties

Warranty information can be obtained by calling AOI Customer service at 949.580.1266.

All orders are shipped FOB Lake Forest.



American Optisurgical Inc.

26902 Vista Terrace
Lake Forest, California 92630

Toll Free: 800.576.1266 (U.S. Only)

Tel: 949.580.1266

Fax: 949.580.1270

email: info@optisurgical.com

website: www.optisurgical.com