



QuickRinse™

The QuickRinse™ Automated Instrument Rinse System standardizes post-surgical care of your reusable instruments and tubing prior to sterilization.



AOI
AMERICAN
OPTISURGICAL
INCORPORATED



The QuickRinse™ Automated Instrument Rinse System — easy and efficient at the same time.

FEATURES

- Designed for rinsing phaco and I/A handpieces, cannulas, vit cutters, reusable tubing and any other micro-surgical instrument with lumen
- Replaces syringe rinsing methods
- Helps reduce damage to instruments from excessive handling
- Automated cycles for consistent rinsing and air drying with a manual override feature
- Compact and lightweight – unit has 8.75” X 8.75” footprint and weighs only 4 pounds
- Bottle, cap, and tubing disconnect easily for sterilization
- **System includes** adapters, nozzles, a cleaning brush and bottle assembly


- **Optional accessories include:**
 - A waterproof foot pedal offers “hands-free” operation
 - Nose cone flush attachment
 - Stellaris™ Irrigation line flush attachment
 - Dual rinse “Y” fitting attachment
 - Filtered air line



System includes:



Optional accessories:



“The QuickRinse is a necessity for our fast-paced surgery center. We move onto other tasks while the QuickRinse is thoroughly flushing the handpieces. I have seen microscopic photos of the insides of phaco tips that were flushed by the QuickRinse and by hand. The difference is astounding. I couldn’t live without it and highly recommend it.”

– A. Busch, Senior Surgical Technician, Overland Park, Kansas



26902 Vista Terrace
Lake Forest, California 92630

Toll Free: 800.576.1266 (U.S. Only)

Tel: 949.580.1266

Fax: 949.580.1270

email: info@optisurgical.com

website: www.optisurgical.com



23511