


# QuickRinse®

## *Instrument Rinse System*



The QuickRinse® from Advanced Optisurgical Inc. is used to rinse Phaco and I/A handpieces, cannulas, vit cutters, reusable tubing and other instruments with lumen. Superior to a disposable syringe, the QuickRinse® provides consistent rinsing pressure and volume regardless of the operator, and eliminates ergonomic issues associated with repetitive syringe use. Using the QuickRinse® also frees up your hands to perform other tasks, greatly improving the speed and efficiency of your reprocessing department.



**AOI**  
ADVANCED  
OPTISURGICAL  
INCORPORATED



# The QuickRinse® Instrument Rinse System automates your rinsing process.

## FEATURES

- **Pushbutton Controls:** Simple user interface consists of only two buttons located on the front panel
- **Multiple Operational Modes:** Timed "auto" cycle (30 seconds) or manual control over cycle time
- **Foot Pedal (Optional):** Offers hands-free operation (convenient for rinsing certain instruments)
- **Quick-Disconnect Fittings:** Enable rapid refilling of the bottle
- **Storage Trays:** Provide convenient location for nylon brush and accessory fittings
- **Colored Tubing:** Aids in proper identification of tubing sets
- **Accessory Fittings and Attachments:** Ensure a compatible connection between the QuickRinse® and your instrument or tube set

### \*NOW AVAILABLE\*

Rinse attachment to fit Centurion and Stellaris handpiece irrigation ports

"The QuickRinse is a necessity for our fast-paced surgery center. We move onto other tasks while the QuickRinse is thoroughly flushing the handpieces. I have seen microscopic photos of the insides of phaco tips that were flushed by the QuickRinse and by hand. The difference is astounding. I couldn't live without it and highly recommend it."

— A. Busch, Senior Surgical Technician, Overland Park, Kansas

## System Includes:

### Accessory Fittings




### Bottle, Bottle Cap, and Tubing



### UL Listed Power Supply



### Spare Quick Disconnect O-Rings



**AOI**  
ADVANCED  
OPTISURGICAL INC.

9 Orchard, Suite 111  
Lake Forest, California 92630

Toll Free: 800.576.1266 (USA only)  
Tel: 949.580.1266  
Fax: 949.215.1199

email: [info@optisurgical.com](mailto:info@optisurgical.com)  
website: [www.optisurgical.com](http://www.optisurgical.com)

This equipment has been verified to comply with the limits for a class B computing device, pursuant to FCC Rules CFR 47 part 15.