

# QuickRinse®

*Instrument Rinse System*

## Accessories and Replacement Parts



425-4003-010



425-4004-010



430-4023-001



750-1100-001



750-2009-001



750-3001-001



750-3008-004



750-3010-001



750-3011-001



750-3013-001



750-3014-001



750-3015-001

If you are located outside of USA, contact your local distributor to purchase any of these accessories or replacement parts.



750-4000-004



750-4001-001



750-5001-001



750-6019-001



750-7000-001



750-7004-001

## Part Number

## Description

425-4003-010	Female Luer Adapters (Qty. 10)
425-4004-010	Male Luer Adapters (Qty. 10)
430-4023-001	Power Cord, (120VAC) USA
750-1100-001	Power Supply, (12VDC) Universal
750-2009-001	Foot Pedal
750-3001-001	Poly Bottle (1000mL)
750-3008-004	Bottle Cap Gaskets (Qty. 4)
750-3010-001	Complete Bottle Assembly (Bottle, Bottle Cap, Air Supply Tubing, Air/Fluid Tubing)
750-3011-001	Bottle Cap with Tubing (Bottle Cap, Air Supply Tubing, Air/Fluid Tubing)
750-3013-001	Air Supply Tubing (Green)
750-3014-001	Air/Fluid Tubing (Blue/White)
750-3015-001	Bottle Cap (with Fittings and Siphon Tube)
750-4000-004	Instrument Interface Tubing (Qty. 4)
750-4001-001	Dual Rinse "Y" Attachment
750-5001-001	Rinse Attachment for Centurion®/Stellaris
750-6019-001	Nylon Brush
750-7000-001	Accessory Fittings Kit [(2) Female, (2) Male Luer Adaptors, (2) 16-Ga. Nozzles, (2) 22-Ga. Nozzles]
750-7004-001	Nose Cone Rinse Attachment



**AOI**  
ADVANCED  
OPTISURGICAL INC.

9 Orchard, Suite 111  
Lake Forest, California 92630

Toll Free: 800.576.1266 (USA only)  
Tel: 949.580.1266  
Fax: 949.215.1199

email: [cs@optisurgical.com](mailto:cs@optisurgical.com)  
website: [www.optisurgical.com](http://www.optisurgical.com)

If you are located outside of USA, contact your local distributor to purchase any of these accessories or replacement parts.